






DX Easy Imaging DR Retrofit Kit - DR Panels Overview

| | | | |
|--|---|---|---|
| Active Area 14" x 17" 358mm X 430mm |  |  |  |
| | PE 1417 W/WL Csi | KM 1417 W/WL Csi | VW 1417 W/WL Gadox/Csi |
| Resolution | 2.304 x 2.800 Pixel | 1.994 x 2.430 Pixel | 2.560 x 3.072 Pixel |
| Pixel Pitch | 150 µm | 175 µm | 140 µm |
| Grayscale | 14 Bit (16.383) | 16 Bit (65.535) | 14 Bit (16.383) |
| Method of exposure | AED (i-Sync2) / DR Trigger | DR Trigger | AED / DR Trigger |
| Limiting Resolution | 3.3 lp/mm | 3.4 lp/mm | 3.5 lp/mm |
| Weight | 3.7 kg (8.16 lbs) | 2.9 kg (6.39 lbs) | 3.1 kg (6.83 lbs) |
| Operating Environment | +5°C to +30°C / 45%-85% rH | +10°C to +30°C / 35%-80% rH | +10°C to +35°C / 30%-85% rH |
| Storage and Transport | -10°C to +40°C / 45%-85% rH | -20°C to +50°C / 20%-90% rH | -15°C to +55°C / 10%-90% rH |

| | | |
|--|--|--|
| Active Area 17" x 17" 430mm X 430mm |  |  |
| | PE 1717 W/WL Csi | KM 1717 W/WL Csi |
| Resolution | 3.072 x 3.072 Pixel | 2.428 x 2.428 Pixel |
| Pixel Pitch | 139 µm | 175 µm |
| Grayscale | 14 Bit (16.383) | 16 Bit (65.535) |
| Method of exposure | AED (i-Sync2) / DR Trigger | DR Trigger |
| Limiting Resolution | 3.4 lp/mm | 3.4 lp/mm |
| Weight | 4.7 kg (10.36 lbs) | 3.6 kg (7.94 lbs) |
| Operating Environment | +5°C to +30°C / 45%-85% rH | +10°C to +30°C / 35%-80% rH |
| Storage and Transport | -10°C to +40°C / 45%-85% rH | -20°C to +50°C / 20%-90% rH |